

Estimation of Ni^{+2}



COMPLEXOMETRICALLY AND GRAVIMETRICALLY

ROHAN SINGH

Chemistry Department, University of Delhi


M.Sc. Chemistry, Group-6, S,No-23

Complexometric Titrations




- ❖ It is a form of volumetric analysis, based on chemical characteristics of complexes and complex forming reactions to estimate relative proportions of metal ions in a sample
- ❖ Appearance or disappearance of coloured complexes dictates the end point
- ❖ EDTA is most widely used as **complexone**
 - Forms highly stable 1:1 complexes due to hexadenticity and chelation
 - Used as disodium salt ($\text{Na}_2\text{H}_2\text{Y}$), other forms are either acidic, basic or insoluble in water
- ❖ Ni^{+2} ions in a solution can be estimated by direct titration against EDTA, using Murexide (or Bromopyrogallol red) as metal ion indicator
- ❖ EBT **cannot** be used as $K_{\text{f}(\text{Ni-EBT})} > K_{\text{f}(\text{Ni-EDTA})}$ (Ni^{+2} is said to block the EBT indicator)

Murexide




- ❖ It is one of the first metal ion indicators, and is **ammonium salt of unstable purpuric acid**. It has 4 protons (1 on each imido group) \rightarrow H_4D^- . Only 2 H can be deprotonated by an alkali hydroxide




❖ Characteristics

- Chromophores \rightarrow Carbonyl ($C=O$) and Imine ($C=N^-$), Auxochrome \rightarrow Imido (batho and hyper chromic shifts)
- Tridentate ligand (2 N and 1 O^-)
- Both **pH** and **pM** sensitive. Hence, pH needs to be maintained using a buffer
- pH range \rightarrow 9-11 (**violet** or **blue-violet** in free state)
- Buffer \rightarrow NH_4OH-NH_4Cl (basic ammonical buffer)
- Complexes \rightarrow Orange (Cu^{+2}), Yellow (Co^{+2} , **Ni^{+2}**), Red (Ca^{+2})
- ❖ Can be used in direct EDTA titration of Ni^{+2} at pH 9-11, end point **yellow** \rightarrow **violet**
- ❖ It has low affinity towards most metal ions, hence used in excess
- ❖ Aqueous solution of murexide are unstable and must be prepared each day




Principle

- ❖ Principle behind complexometric estimation of Ni^{+2} is that it forms stable complex with EDTA, which is more stable than its complex with indicator (here murexide) ($K_{f(\text{Ni-EDTA})} > K_{f(\text{Ni-In})}$)
- ❖ If to an unknown sample of M^{+n} , suitable buffer (to maintain pH) and metal-ion indicator were added, indicator being a ligand forms complex with metal ions imparting colour to the solution (stability and colour depends on M^{+n} and pH)
- ❖ It can then be titrated against standard EDTA, which forms $[\text{M-EDTA}]^{n-4}$ complex owing to its high stability. Near end point, EDTA progressively displaces M^{+n} ions from $[\text{M-In}]$ complex
- ❖ At the end point, all of the M^{+n} ions were complexed by EDTA, leaving free indicator in the solution which causes a colour change



- ❖ H^+ ions produced are neutralized by buffer



4. Complexometric Reactions

Requirements

❖ Chemicals Required

- Standard -> ZnO as primary standard
- Complexone -> Disodium EDTA
- Buffer -> NH_4Cl + NH_4OH solution (ammonical buffer)
- Indicators -> EBT and Murexide

❖ Apparatus Required

- Conical flasks and beakers
- Volumetric flask and measuring cylinder
- Test tubes
- Pipette
- Burette
- Burette stand and white tile


Procedure

1. Preparation of standard ZnO solution (100 mL)

Calculated amount of ZnO (s) + Few drops of HCl (to completely dissolve ZnO in water) + Distilled water (up to 100 mL mark) --> Mixture was shaken well to get homogenous solution




2. Standardization of EDTA solution



3. Estimation of Ni^{+2}

- Indicator \rightarrow 0.1 g murexide + 10 g KNO_3 , grind well and use 50 mg of mixture in each titration
- Pipette 25 mL of given Ni^{+2} solution in a conical flask and dilute to 100 mL (with deionised water)
- Add 10 mL of 1 M NH_4Cl + conc NH_4OH buffer followed by 50 mg of indicator mixture
- If solution turns orange, add more conc NH_4OH until solution turns clear yellow (pH~10)
- Titrate against standard EDTA solution dropwise, until colour of solution changes from **yellow** \rightarrow **violet** at the end point (If at any intermediate stage, solution becomes orange, conc NH_4OH must be added to make it clear yellow)
- Ni^{+2} complexes slowly with EDTA, hence it must be added **dropwise near the end point**



Calculations


Let normality of standard ZnO solution prepared be N_{ZnO} , then

❖ Normality of given EDTA solution



- $N_{\text{EDTA}} \times V_{\text{EDTA}} = N_{\text{ZnO}} \times V_{\text{ZnO}}$ ($V_{\text{ZnO}} = 10 \text{ mL}$)
- $N_{\text{EDTA}} = N_{\text{ZnO}} \times (10 / V_{\text{EDTA}})$ (V_{EDTA} obtained from step 2)


❖ Normality of Ni^{+2} in given solution



- $N_{\text{Ni}^{+2}} \times V_{\text{Ni}^{+2}} = N_{\text{EDTA}} \times V_{\text{EDTA}}$ ($V_{\text{Ni}^{+2}} = 25 \text{ mL}$)
- $N_{\text{Ni}^{+2}} = N_{\text{EDTA}} \times (V_{\text{EDTA}} / 25)$ (V_{EDTA} obtained from step 3)

Gravimetry

- ❖ It's a form of quantitative analysis by weight, and involves isolation and weighing of specie being determined (or compound of that specie)
- ❖ **Principle** -> If an element or radical under estimation forms a pure and stable compound with definite composition that can be readily isolated and weighed, then the amount of that element or radical in the test sample can be calculated just by relative atomic masses of constituent elements



- ❖ Only significant disadvantage of this analysis technique is it's time consuming nature, which limits it's use in many places

Scope

❖ Advantages

- **Absolute method**, involving only direct measurements without any form of calibration
- Highly accurate and precise, with modern weighing balances
- Regular checking of possible error sources
 - Filtrate -> can be checked for completeness of precipitation
 - Precipitate -> can be checked for impurities
- **Mostly uses inexpensive apparatus** (the most costly being muffle furnace and platinum crucible in rare cases)



❖ Applications

- Analysis of standards for testing and calibration of instrumental techniques
- Analysis which require high accuracy
- Estimation of chemical composition, purity and thermal stability of compounds



Precipitation

It is the most widely used isolation method in gravimetric analysis. The element (or radical) under analysis is transformed into an insoluble precipitate, which can then be isolated

IDEAL PRECIPITATE

- ❖ It should be readily convertible to a pure compound with definite chemical composition, via simple chemical processes like
 - Ignition
 - evaporation with a suitable liquid, etc
- ❖ It should be highly insoluble, which on practical grounds means that
 - Very low K_{sp} (\ll Ionic Product)
 - Amount of it left in solution should be less than least count of weighing balance used
- ❖ Physical nature should be such that it can be
 - Readily separable by filtration (large particle size)
 - Washed free of soluble impurities (particle size unaffected or at least not diminished by washing)

Precipitate Contamination

It's not necessary that precipitate that separates out of solution is always pure. Type and degree of contamination depends upon it's nature and condition of precipitation

CO-PRECIPIATION

Contamination of ppt by substances which are normally soluble in mother liquor (co -> together)

❖ Surface Adsorption

- Adsorption of contaminants over surface of ppt particles exposed to solution
- Greatest for gelatinous (large surface area) ppt, minimum for macro-crystalline (large particle size) ppt

❖ Internal Contamination

- Mostly occurs during fast building of precipitate from primary particles (coalescence)
- Occlusion -> physical entrapment of contaminants adsorbed over primary particles (during fast coalescence) within ppt
- Inclusion -> same as occlusion, but contaminants now occupies lattice sites within the ppt

POST-PRECIPIATION

- Precipitation that occurs over the surface of main precipitate *after (or post)* it's formation
- Increases with time main ppt is left in contact with mother liquor
- Ex $\text{Mg}(\text{C}_2\text{O}_4)$ post precipitates over $\text{Ca}(\text{C}_2\text{O}_4)$ ppt

Digestion

Digestion is carried out by allowing ppt to stand for 12-24 hours or warming ppt for some time in contact with mother liquor

PURPOSE


- reduce contamination
- increase particle size and complete precipitation

EFFECT

- ❖ Small particles
 - have high solubility and tend to pass in solution and re-deposit over larger particles
 - Co-precipitation over minute particles is hence eliminated
- ❖ Rapidly formed crystals
 - have irregular shape and high surface area, which on digestion becomes regular and dense
 - Result is decrease in surface area and adsorption



RESULT


- Decrease in overall co-precipitation
- Increase in particle size and ease of filtration



DMG

- ❖ Dimethylglyoxime { $\text{H}_3\text{C} \cdot \text{C}(\text{:NOH}) \cdot \text{C}(\text{:NOH}) \cdot \text{CH}_3$ } was discovered by L.Tschugaeff, and applied by O.Brunck in estimation of Ni^{+2} in steel
- ❖ Has 2 ionisable H $\rightarrow \text{H}_2\text{DMG}$
- ❖ Very slightly soluble in water (0.40g/L) \rightarrow used as 1% solution in 90% ethanol (rectified spirit) or absolute ethanol
- ❖ Can be used as disodium salt ($\text{Na}_2\text{C}_4\text{H}_6\text{O}_2\text{N}_2 \cdot 8\text{H}_2\text{O}$), which is soluble in water
- ❖ Bidentate ligand, and forms a bright red $[\text{Ni}(\text{HDMG})_2]$ complex with Ni^{+2}
- ❖ Large excess of DMG should be avoided because
 - It can itself precipitate due to low solubility in water
 - $[\text{Ni}(\text{HDMG})_2]$ ppt can dissolve in high ethanol content (over 50%)
- ❖ Interference from Fe^{+3} , Al^{+3} , Cr^{+3} and Bi^{+2} is prevented by adding soluble tartrate or citrate, with which they form soluble complexes






- ❖ Ni^{+2} solution is made faintly acidic (by dil HCl), to prevent it's hydrolysis (as $\text{Ni}(\text{OH})_2$ is very insoluble, $K_{\text{sp}} \sim 5.85\text{E}-16$) and heated ($70-80^\circ\text{C}$) to increase complexation rate
- ❖ $[\text{Ni}(\text{HDMG})_2]$ ppt is **soluble** in
 - Free mineral acid (even as low as produces by reaction)
 - Solutions having > 50% ethanol (by volume)
 - In concentrated ammonia
- ❖ Therefore...
 - Large excess of DMG containing ethanol is avoided
 - **slight** excess of dil NH_4OH or ($\text{CH}_3\text{COONH}_4 + \text{CH}_3\text{COOH}$) buffer is added to neutralise free mineral acids
 - CH_3COONa in place of $\text{CH}_3\text{COONH}_4$ should be used in presence of Mn, Co and Zn

15. Precipitation Reaction

Requirements



DESICCATOR

❖ Chemicals Required

- Ni^{+2} solution
- Dilute HCl solution (1:1)
- DMG reagent (1% solution in 90% or absolute ethanol)
- Dilute NH_4OH solution

❖ Apparatus Required

- Crucible
- Desiccator (sealable container with **desiccants** which absorb moisture from air, ex CaCl_2 , CaSO_4 . Used to sustain state of dryness)
- Weighing balance
- Beakers (100 and 500mL)
- Test tubes
- Measuring cylinder
- Glass rod



BALANCE



CRUCIBLE

Procedure

❖ Preparation of Ni^{+2} solution (200 mL)

0.3-0.4g $(\text{NH}_4)_2\text{SO}_4 \cdot \text{NiSO}_4 \cdot 6\text{H}_2\text{O}$ (Ammonium Nickel Sulphate) in 500 mL beaker + 5 mL dil HCl
(1:1, to prevent hydrolysis) + **dilute to 200 mL with distilled water**

❖ Precipitation of $[\text{Ni}(\text{HDMG})_2]$

- **Heat Ni^{+2} solution to 70-80°C, and add slight excess of DMG reagent** (at least 5 mL for every 10 mg Ni^{+2})
- **Immediately add dil NH_4OH solution dropwise** (directly to solution rather down the sides of beaker) **with continuous stirring until precipitation occurs, and finally in slight excess**
- **Allow to stand on steam bath for 20-30 minutes, and then check for complete precipitation**
- **Allow ppt to stand for an hour, cooling at the same time**
- **Filter cold solution through a porcelain filtering crucible** (preheated to 110-120°C)
- **Cool the Nickel Dimethylglyoximate (bright red) precipitate in desiccator**

Drying and Weighing

- Wash $[\text{Ni}(\text{HDMG})_2]$ ppt with cold water (until free from soluble salts like chlorides)
- Dry at $110\text{--}120^\circ\text{C}$ (or 150°C for high accuracy) for 45-50 minutes. This volatises any reagent that may have carried down with ppt
- Cool in desiccator (moisture free environment) and weigh
- Repeat the drying until constant weight is achieved



Calculations

❖ % of Ni in $[\text{Ni}(\text{HDMG})_2]$

- Molecular Formula $\rightarrow \text{Ni}(\text{C}_4\text{H}_7\text{O}_2\text{N}_2)_2$
- Relative atomic masses of constituents (u) $\rightarrow \text{Ni} = 58.69, \text{C} = 12, \text{H} = 1, \text{O} = 16, \text{N} = 14$
- Total Molecular Mass (u) $\rightarrow 288.69$
- % of Ni $\rightarrow (58.69\text{u} / 288.69\text{u}) * 100 = 20.32\%$

❖ Mass of Ni^{+2} in test sample ($M_{\text{Ni}+2}$)

Let mass of $[\text{Ni}(\text{HDMG})_2]$ measured experimentally (constant value) $\rightarrow M_{\text{Ni-DMG}}$

- Then, $M_{\text{Ni}+2} = \text{Mass of } \text{Ni}^{+2} \text{ in } [\text{Ni}(\text{HDMG})_2]$
- $M_{\text{Ni}+2} = (\% \text{ of Ni in complex} / 100) * M_{\text{Ni-DMG}}$
- $M_{\text{Ni}+2} = (20.32 / 100) * M_{\text{Ni-DMG}}$
- $M_{\text{Ni}+2} = 0.2032 * M_{\text{Ni-DMG}}$

Comparison

Complexometry	Gravimetry
Reasonably accurate (depends upon stability of M-In and M-EDTA complex)	Highly accurate (depends upon nature of specie being isolated)
Quick estimation	Consumes lot of time
Relative method (involves titration against standards)	Absolute method (no calibration required)
Highly error prone (pH sensitivity, low stability of complexes, indicator error etc)	Error sources (filtrate and precipitate) can be easily checked
Can produce a lot of waste	Waste is mostly limited to reagents and isolated compound
Used to estimate relative proportions of metal ions in the solution	Used in high accuracy demanding areas, like analysis of standards

References

- ❖ Vogel's textbook of quantitative chemical analysis
- ❖ <https://en.wikipedia.org>
- ❖ <https://chemistry.stackexchange.com>
- ❖ <https://www.chemeuropa.com/en/encyclopedia>